

## • ABOUT US •

Driven by the passion for innovation, we at Dr Trust endeavour to provide our customers with the latest medical inventions with an objective to promote good health and wellness all around the world. All the medical devices and health monitors provided by Dr Trust are supported by accurate, latest and ground breaking technologies, innovated at our headquarter in NY, USA. All our products adhere to the most stringent CE and FDA guidelines and are strongly recommended by doctors and health practitioners. Our products are designed in the utmost exemplary ways to ensure that their accuracy and convenience are unrivalled. The ease of their use and operation makes them even more suitable for users of all age groups.

Dr Trust strives to enhance the quality of lifestyle by providing with the most trusted and innovative health care and wellness products. Being a renowned global leader in health care products, Dr Trust ensures that our technically efficient team works dynamically and tirelessly to provide the best of the medical devices to our clients. The products that we have to offer are suitably designed for use at homes, laboratories and hospitals.

Our ground breaking solutions allow you to monitor your health in the easiest ways possible. In today's era when all of our lives are too hassled to handle, it becomes a bit difficult to pay attention to our health. But it has now become easier with the coming of the monitoring devices which can be conveniently used at homes and even on the go.

**We bring to you a variety of best self medical devices, trusted and used by Doctors, medical professionals and home users all over the world.**

## Dr Trust

### HomeSpa Facial Steamer & Vaporizer – 901 (Pink)

## • QUICK STARTUP GUIDE •

### IMPORTANT

Arroma oil is not included in the package.

#### Step 1

Remove the steam mask that covers the chamber.

#### Step 2

Add about 40 ml water into the chamber with measuring cup.

#### Step 3

Put back the steam mask over the slot and turn it on.

#### Step 4

Turn the knob to 1-2 position and let it heat up.

#### Step 5

Enjoy evaporated steam and turn the switch off. once you are finished.

**Note:** You can add few drops of aromatic oil on a felt strip or pad then place it inside the aroma tray or the steam nozzle cap. You can also add any herbs of your choice in it.



## • INTRODUCTION •

Dr Trust Home Spa facial steamer & vaporizer is simply designed to emit steam that specifically targeted every part of your face for a proper treatment like spa at home. It is a stylish and exquisite craftsmanship with a power supply cord. It comes with two cones - the wide one fits entirety of your face for proper steaming and the narrower one emits targeted steam for your nose. It is a high-tech steamer supposed to emit regular steam to moist the skin and renew cells to make skin fleshy and elastic.

## • RIGHT USAGE •

### How to put the facial steamer to work?

1. Fill the Water Tank  
The steamer comes with a water tank. Fill the tank with lukewarm distilled water instead of regular tap water.
2. Turn on the Steamer  
Turn on the steamer and allow the water to heat up so that steam starts coming out from the facial steamer.
3. Steam Your Face  
Allow the steam to fall over your face from a certain distance. You can control the steam delivery as per your liking by adjusting the control knob.
4. Steam Timing  
Do this for about 3 to 10 minutes, depending on your skin's type, needs and reaction.
5. Use a Face Mask  
Apply a face mask to reap all the benefits of your steam session.
6. Moisturize Your Skin  
With any skin care routine, finish with a hydrating moisturizer.

## Correct steam timings for various skin types

Its steam lasts about 8 minutes. However, keep in mind that you can receive many of the same benefits if you stop after 5 minutes. Do not steam your face for much longer than 5 minutes, especially if you have acne or other skin issues. Steaming causes the face to swell, and it can exacerbate acne if done for too much time.

### How often should you use?

Listen to your skin to determine the frequency that is best for you. Generally, experts recommend steaming your face once a week. The frequency with which you should steam your face will be dependent on your skin, its condition, and level of sensitivity.

## • TIPS •

### Using the Facial Steamer Safely & Effectively

- Before beginning, follow daily face skin care routine to clean the skin properly.
- Lightly pat your favorite moisturizer or cream around the eye area.
- If more steam is desired, place a towel over your head that so it covers the mask and sides of the steamer.
- Close your eyes and relax as the soothing mist begins deep cleansing action.
- Remove your face occasionally and blot with a cool, damp towel.
- Rinse your face with cool water immediately after steaming.



• **CAUTION** •

1. If you have heart expansion disease, serious rosacea, eczema, or psoriasis, please do not use this product.
2. Pregnant women and patients with serious diseases, such as in the use of the products need to add herbs and essential oil before, please consult a physician.
3. Once appear the skin inflammation or discomfort, immediately stop using essential oils and herbs.
4. Often clean and detergent can ensure the best effect and extend the product life
5. Always unplug the steamer immediately after use.
6. Do not place in or drop into water or other liquids. Never use while bathing.
7. Do not reach for an appliance that has fallen into water or other liquids. Firstly, unplug it immediately.
8. This appliance is not intended for use by persons (including children) with reduced physical, sensory or mental capabilities, or lack of experience and knowledge, unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
9. Children should be supervised to ensure that they do not play with the appliance.
10. Do not use this appliance near bathtubs, showers, basins, or other vessels containing water.

• **SPECIFICATIONS** •

- |                                |  |
|--------------------------------|--|
| 1. 110-240V 50/60Hz 100w power | 4. 0/1/2 Switch setting with LED indicator light |
| 2. PTC heating                 | 5. Overheating protection prevents the dry burn  |
| 3. CE / ROHS certification     | 6. Flexible as it has an adjustable steam size   |

• **MAINTENENCE** •

1. After each use, it should be cleaned properly to avoid oil and other dirt deposition.
2. Always turn off and unplug it from the electrical outlet before you start cleaning.
3. Let it cool down, then empty excess water, and wipe dry interior of the facial mask and rim around the opening.

• **CUSTOMER SUPPORT** •

**CONTACT ADDRESS**

**USA**

Nureca INC.USA  
276 5th Avenue, Suite 704-397,  
New York (NY) - 10001, USA

**INDIA**

**Corporate Office (Mumbai)**

Nureca Limited  
128 Gala Number Udyog Bhavan,  
1st Floor Sonawala Lane, Goregaon East  
Mumbai City Maharashtra 400063

**Contact us**

India: +91-7527013265 /+91-9356658436

Website: [www.drtrust.in](http://www.drtrust.in)

Corp Website: [www.nureca.com](http://www.nureca.com)

Email: [customercare@nureca.com](mailto:customercare@nureca.com)

Connect with us on social networks

Facebook: @drtrust

Instagram: @drtrustisin

Youtube: NurecaUsa

COPYRIGHT©2021 NURECA LTD ALL RIGHTS RESERVED

